

# HF RADDAM Status

(Radiation Damage Monitoring System)

Yalcin Guler<sup>1</sup>, Jean-Pierre Merlo<sup>2</sup>, Taylan Yetkin<sup>2</sup>


1-University of Cukurova

2-University of Iowa

# Environment for Fibers



- HF detectors will experience very high radiation levels (100 Mrad at  $\eta=5$  with  $500 \text{ fb}^{-1}$  )
- After  $\sim 10$  years of operation the fibers will lose their light attenuation by half.
- The degradation will occur gradually and a monitoring tool is necessary.

# RADDAM principle [3]




- Inject light through a capillary tube in a 2.5 m long fiber (two ends polished). For HF this length corresponds to a regular fiber length read by PMT.
- Reflection occurs at the two ends, signal S1 at the entrance and S2 at the far end. S2 coming 25 ns later.
- The ratio  $R = S1 / S2$  is related to the fiber transparency, and its value depends on the accumulated dose and post irradiation time  $t$  for recovery.
- The fibers in HF are grouped in towers. Raddam fibers are equipped in 7 eta towers of 8 HF wedges.
- HF detected light (fibers+ PMT) has a spectrum centred at 450 nm

# Radiation Damage Monitoring Fiber





# HF Raddam device





The device use a pulsed nitrogen laser (337 nm) and it is wavelength shifted into the blue region ( 440 nm) by a 2 cm long scintillating fiber positioned in a feed-through connector. This light is then distributed into fibers in four wedges per HF. Within a wedge, 7 towers (spanning all rapidity regions) are instrumented with these sample fibers. Each fiber is readout by the PMT that serves that particular tower.[6]

# Radiation damage in irradiated quartz fibre [1][2]



$\alpha$  and  $\beta$  parameters versus  $\lambda$  for qq and qp fibres corresponding to the fit of irradiation data up to 1.25 Grad.

# Damage Recovery parameters [1][2]



Recovery parameters  $\gamma$  and  $\eta$  versus  $\lambda$  for qq and qp fibres versus wavelength, after 205 and 605 h of irradiation.

# Damage and recovery of irradiated quartz fibres [1][2]

Radiation damage (Decrease of signal)

$$I(\lambda, D)/I(\lambda, 0) = \exp[-A(\lambda, D).L/4.343] \rightarrow S1/S2$$

Optical attenuation is defined as

$$A(\lambda, D) = \alpha(\lambda) [D/D_s]^\beta$$

where

D: Accumulated dose

D<sub>0</sub>: Reference dose

$\alpha(\lambda)$  and  $\beta(\lambda)$ : Characteristic parameters which define radiation hardness of a given fiber (must be found empirically).

A in dB/m, D in Mrad, L in m

Recovery (Increase of transm. signal)

$$\frac{I_D(\lambda, t)}{I(\lambda, D)} = \exp \left[ A(\lambda, D) (L/4.343) \left( \frac{\gamma t^\eta}{t_{irr}^\eta + \gamma t^\eta} \right) \right]$$



# Estimate of radiation damage and recovery versus dose in HF quartz fibers

$t_{\text{irr}}$ (days)	1	2	3	5	10	20	50	100	150	200	2000
<b>Dose</b> (Mrad)	<b>0.5</b>	<b>1.0</b>	<b>1.5</b>	<b>2.5</b>	<b>5.0</b>	<b>10</b>	<b>25</b>	<b>50</b>	<b>75</b>	<b>100</b>	<b>1000</b>
<b>A</b> (dB/m)	<b>0.156</b>	<b>0.211</b>	<b>0.252</b>	<b>0.316</b>	<b>0.428</b>	<b>0.581</b>	<b>0.869</b>	<b>1.179</b>	<b>1.410</b>	<b>1.600</b>	<b>2.754</b>
<b>Signal loss (%)</b>	<b>5.7</b>	<b>7.7</b>	<b>9.1</b>	<b>11.3</b>	<b>15.0</b>	<b>19.8</b>	<b>28.0</b>	<b>36.1</b>	<b>41.5</b>	<b>45.6</b>	<b>65.0</b>
<b>Increase aft.1h (%)</b>	<b>1.7</b>	<b>1.9</b>	<b>2.1</b>	<b>2.3</b>	<b>2.7</b>	<b>3.0</b>	<b>3.6</b>	<b>4.0</b>	<b>4.3</b>	<b>4.5</b>	<b>4.1</b>
<b>Increase aft.8h (%)</b>	2.5	3.0	3.3	3.8	4.4	5.1	6.2	7.1	7.7	8.1	7.4
<b>Increase aft.1d (%)</b>	3.0	3.7	4.1	4.7	5.6	6.6	8.1	9.0	9.8	10.9	10.2
<b>Increase af.10d(%)</b>	<b>4.0</b>	<b>5.1</b>	<b>5.8</b>	<b>6.8</b>	<b>8.5</b>	<b>10.4</b>	<b>13.4</b>	<b>16.1</b>	<b>17.9</b>	<b>19.2</b>	<b>19.4</b>





This table shows the increase of 450 nm signal at different dose and post d.t. time t (dose rate 0.5 Mrad/day)

# Irradiation test history [3]

PS Irradiation at high dose rate 2.6 Mrad/hour, 250 Mrad in 4 days



# Event Selection



**radTS** : time slice id where S1 measured. radTS = 3 for HF(+), 4 for HF(-)

**2Q** : sum of the fC in radTs and radTS+1.


**4Q** : sum of the fC in radTs-1, radTS, radTS+1, and radTS+2.

**Selection criteria:**


$2Q/4Q > 0.95$  &  $Q[\text{radTS}] > 50 \text{ fC}$  &  $Q[\text{radTS}+1] > 50 \text{ fC}$   
 (50 fC requirement is used after pedestal subtraction)

On average ~20-30% of the total events remain.

Ch 51 HF(40,71,1)




Ch 51 HF(40,71,1)




After event selection S1 and S2 distribution are fitted with gaussian fit.


Ch 51 HF(40,71,1)



Ch 51 HF(40,71,1)



Ch 51 HF(40,71,1)




We considered the following error sources:

- 2% QIE quantization error (i.e. effect of binning)[4]
- 2.5% QIE-to-QIE spread [5].

The resulted error estimation is 3.2%.

### Ch 51 HF(40,71,1)




2010 Runs

$$\langle S1(t) \rangle / \langle S1(0) \rangle$$

Monitoring of laser light intensity

### Ch 51 HF(40,71,1)




2010 Runs

$$\langle S2(t) \rangle / \langle S2(0) \rangle$$

Laser light intensity + Raddam effect

# Ch 51 HF(40,71,1)





$$R(t)/R(0) = [(\langle S1(t) \rangle / \langle S2(t) \rangle) / (\langle S1(0) \rangle / \langle S2(0) \rangle)]$$

Normalized Raddam effect

# DQM Update

Plan to update on DetDiag Raddam offline analyzer :

- Implement for each channel and Run S1, S2 distributions and S1/S2
- Update plot view for short-term problems lasers, PMT, QIE and timing.
- Provide a stand-alone ROOT macro for trend analysis based on DQM outputs.



# WHAT COULD BE IMPROVED

- A good laser with more precise timing and less aging presently we select only 20-30 % of data
- Install spares of long fibers going from laser shifter to HF±
- Replace SMA connectors by ST ones at the level of laser bench (radial positioning of fibers)
- Improved event selection.
- Produced one example DQM plot, will work on updating DetDiag for RADDAM runs.
- Final goal is to provide a simple ROOT macro for trend analysis based on DQM outputs.

# Conclusion

- Raddam device is working in 55 channels /56
- We have reliable analysis of Raddam data to be installed in DQM
- Can be used to monitor:
  - Radiation damage of fibers,
  - Laser intensity
  - Stability of electronics, quality of timing in 55 channels

# References

[1] I.Dumanoglu et al. NIM A 490 (2002) 444-455

[2] K. Cankocak et al., CMS note 2007/03 submitted to NIM A. 585(2008) 20-27

[3] <http://indico.cern.ch/getFile.py/access?contribId=7&resId=1&materialId=slides&confId=38625>

[4] <http://lss.fnal.gov/archive/2010/pub/fermilab-pub-10-210-ppd.pdf>

[5] <https://twiki.cern.ch/twiki/bin/viewauth/CMS/ChargeInjector>

[6] CMSNOTE 2006-044