

Module testing status and X-ray tests

Suvadeep Bose

University of Nebraska Lincoln

(On behalf of the SIDET team)

- Status at SIDET
- Status of test module test stands
- Person power at SIDET and future plans
- Work with the new TBM08 on HDI
- X-ray test stands – at UIC and KU

- ☐ We have set up the two cold boxes from UIC
- ☐ There are two PC's (Ubuntu linux) and psi46 softwares installed on them
- ☐ The chillers are set up
- ☐ Cold boxes are connected to the N₂ source
- ☐ New HV (Keithley) are in place
- ☐ *El Commandante* is talking to the JUMO controller in the Cold boxes
- ☐ We have succeeded to talk to the Cold boxes and to the DTB at the same time
- ☐ We have two test stands set up – full commissioning yet to happen

☐ News on chips:

- We received one bare 1x1 BPix sensor last month which we tested and sent to Purdue.
- Currently tests are being performed on
 - one 1x1 bare single ROC
 - one 1x1 bare single ROC for FPix bonded to a sensor

Test station 1




Test station 2




What is available at SIDET ?

1 HDI with TBM08





1 Adapter board



1 Prototype Module Carrier board



(We have 5 of them in US, materials have been purchased for 50 and they are being machined at UIC)



In addition we have


- 4 Digital Test Boards (DTB)
- 12 Single ROC adapter cards
- 2 x 68-pin ribbon cables (standard SCSI cables can be used)

- We tried placing the prototype carrier board on the cold plate of the Cold box
- The carrier board fits well but there is a **thermometer** at the edge which will prevent the fourth carrier board to sit properly on the plate.




- We have to make arrangements to move it to somewhere else

- We also tried placing the HDI inside the prototype carrier board
- It fits fine.




- Person power at SIDET:
 - Douglas Berry (Postdoc, UIC)
 - Suvadeep Bose (Postdoc, UNL)
 - Rafael Coelho Lopes De Sa (Postdoc, Fermilab)
 - Cecilia Gerber (UIC)
 - Daniel Sandoval (PhD student, UIC)
 - John Stupak (Postdoc, Purdue Calumet)
 - Paul Turner (PhD student, UIC)
 - Marco Verzocchi (Fermilab)
- We meet weekly on Mondays at 3 PM (22:00 GVA) in Wilson Hall to discuss the status of work
- In next 2/3 weeks we plan to fully commission the cold boxes, calibrate them and determine the operating conditions.
- In longer term we plan to understand the psi46 software and get involved in its development

- ❑ We wanted to test new firmware that Robert (KU) developed. We have a new Token Bit Manager (TBM08) wire bonded to an HDI without a sensor.



- ❑ We have observed the same problem that was seen by others while testing the TBM.

We tried sending a signal and tried to look at the readout on an oscilloscope. As the test was done with an HDI without sensors the readout will be empty.



- ❑ We were trying to readout the TBM with the **psi46test** with the new firmware (1.13)
- ❑ We tried to **decode** the idle pattern with the DESERializer sending a *single trigger*.

When TBM generates an **idle pattern** like 0000111100001111

Running DESER we saw that:

dread 0

#samples: 0 remaining: 0

(The dread 0 shows the decoded channel 0)

> dread 1 ((shows a test pattern)

#samples: 10 remaining: 99988

C155=1100000101010101

C156=1100000101010110

C157=1100000101010111

C158=1100000101011000

C159=1100000101011001


C15A=1100000101011010

C15B=1100000101011011

C15C=1100000101011100

C15D=1100000101011101

C15E=1100000101011110



If the DESER receives the idle pattern it does not record anything, thus the 0 samples.
When a trigger comes - the data is recorded.

❑ In first look – the DESER seems to be working correctly but further tests will be done at UNL and Fermilab.

- ❑ 0-50kV, 0-1mA X-Ray tube mounted on rotary mount for both direct and indirect X-Ray beam
- ❑ Fluorescent foils (Cu, Mo, Ag, In, Sn, Tb) and custom computer controlled jig for target selection
- ❑ Thermoelectric cold plate with temperature controller (air cooled) w/ support for single module or single ROCs. Capable of chilling down to $\sim 5^{\circ}\text{C}$
- ❑ Software integration into elComandante to allow “push button” testing



- ❑ 0-50kV, 0-1mA X-Ray tube mounted on rotary mount for both direct and indirect X-Ray beam
- ❑ Fluorescent foils (Cu, Mo, Ag, In, Sn, Tb)
- ❑ Thermoelectric cold plate with temperature controller (into water bath with radiator).

